

# HYDROMERSE®

## LED LIGHT MASK

### HYDROMERSE LED LIGHT MASK

Hydromerse LED light therapy face, neck & hands mask can be used by professionals and clients at home to:

- Improve skin tone and texture
- Smooth fine lines and wrinkles
- Increase collagen production
- Reduce acne and prevent breakouts

### Instructions For Use

1. Ensure face and neck are clean and dry, we recommend using Hydromerse Cosmeceuticals Dragons Blood facial cleanser.
2. For maximum benefit apply Hydromerse Hydrating recovery mask – following instructions for use
3. Connect your Hydromerse LED light mask to the remote control – you may also wish to connect the neck LED mask at the same time.
4. Plug the black power adapter into the remote-control unit and connect to power source
5. The Unit is set to 30 minutes when switched on. You can change this using the remote control by clicking the down or up arrow labelled “Time” (3) – maximum 60 minutes.
6. Select the colour of the Hydromerse LED mask you wish to use by pressing the “Colour/Start” button (4).



Blue

Cyan

White



Purple

Red

Green

Yellow

C0 – Off

C1 – RED

C2 – GREEN

C3 – BLUE

C4 – YELLOW

C5 – CYAN

C6 – PURPLE

C7 – WHITE – All colours

C8 – AUTO CHANGE

7. Adjust the brightness by clicking the “Brightness” button display reads A1-A9 intensity. (5)
8. Select “Flicker” button to choose if you’d like the LED lights to flicker display reads F1-F4. (2)
9. Select “Microelectronics” display reads B0 – B9 (B0 is off) B1-B9 is the strength of the micro electric pulses. (1)
10. Find a suitable, comfortable position – lying down is best, and fit the Hydromerse LED LIGHT Mask over your face using the elasticated straps.
11. Once your time has elapsed remove the mask and clean prior to next use if the mask has product on it.
12. Remove your Hydromerse Hydrating Recovery Mask and apply Hydromerse Cosmeceuticals as desired.



Detachable neck mask may be used for the hands or anywhere on the skin

Repeat each day for continued benefits 😊



**Microcurrent Contact**  
Can strengthen the absorption of skin nutrition



**Silicone Pads**  
Intimate eye protection



**192\* LED Lamp Beads**  
Produce more rapid and deep beauty effect



**Neck Connector**  
Can be used separately or together

### LED light Colour and Benefits



**Red light: To treat fine lines & wrinkles (630nm)**

Red light stimulates skin cells to produce collagen. These cells are responsible for skin firmness and elasticity. It improves the skin texture, shrinks pores and reduces fine lines and wrinkles, resulting in a smooth, youthful skin.



**Green light: To improve tone & pigmentation (520nm)**

Green light limits the excess of melanin production in the lower lever of the dermis, preventing it to reach the surface of the skin. It breaks down existing melanin clusters that result in the reduction of age-spots and hyperpigmentation, leaving a smooth and even skin tone.



**Blue light: To control acne & blemishes (470nm)**

The wavelength of blue light concentrates on acne-causing bacteria. This light reaches Porphyrins both above and below the skin, these organisms create breakouts and inflammation on the skin. This process results in the reduction of acne, and calming inflammation and redness, leaving the skin clearer, brighter, and smoother.



**Yellow light:** To improve clarity & texture (590nm)

This wavelength stimulates circulation and improves the cells ability to eliminate toxins by increasing the exchange of oxygen within skin cells. This process improves tone, clarity, and texture, resulting in a glowing complexion.



**Purple light:** To treat redness (415-620nm)

This light combines the benefits of blue and red light. It has an antibacterial effect that helps treat acne, while reducing fine lines and wrinkles as it stimulates the production of collagen.



**Cyan light:** To reduce scarring & inflammation (490nm)

This wavelength assists in the reduction of swollen capillaries. It promotes healing and soothes inflamed and damaged skin, improving clarity and skin tone.



**White light (all LED's):** For the deepest skin penetration and anti-ageing

This wavelength penetrates deep into the skin cells, accelerating the metabolism of active tissue, improving the production of collagen and elastin.



Website for further details: [www.hydromerse.com](http://www.hydromerse.com)

Instagram/Facebook/Twitter - #hydromerse @hydromerse  
LinkedIn: [Hydromerse](https://www.linkedin.com/company/hydromerse) Pty Ltd